

兔抗 DEFB1 多克隆抗体

中文名称：兔抗 DEFB1 多克隆抗体

英文名称：Anti-DEFB1 rabbit polyclonal antibody

别名：BD1; HBD1; DEFB-1; DEFB101

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DEFB1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Defensins form a family of microbicidal and cytotoxic peptides made by neutrophils. Members of the defensin family are highly similar in protein sequence. This gene encodes defensin, beta 1, an antimicrobial peptide implicated in the resistance of epithelial surfaces to microbial colonization. This gene maps in close proximity to defensin family member, defensin, alpha 1 and has been implicated in the pathogenesis of cystic fibrosis.
Applications:	ELISA, IHC
Name of antibody:	DEFB1
Immunogen:	Synthetic peptide of human DEFB1
Full name:	defensin beta 1

Synonyms:	BD1; HBD1; DEFB-1; DEFB101
SwissProt:	P60022
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-100

