

兔抗 DEFB4A 多克隆抗体

中文名称: 兔抗 DEFB4A 多克隆抗体

英文名称: Anti-DEFB4A rabbit polyclonal antibody

别名: BD-2; SAP1; DEFB2; DEFB4; HBD-2; DEFB-2; DEFB102

抗原: DEFB4A

储存: 冷冻 (-20°C)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Defensins form a family of microbicidal and cytotoxic peptides made by neutrophils. Members of the defensin family are highly similar in protein sequence. This gene encodes defensin, beta 4, an antibiotic peptide which is locally regulated by inflammation.
Applications:	ELISA, IHC
Name of antibody:	DEFB4A
Immunogen:	Synthetic peptide of human DEFB4A
Full name:	defensin, beta 4A

Synonyms :	BD-2; SAP1; DEFB2; DEFB4; HBD-2; DEFB-2; DEFB102
SwissProt:	O15263
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

