

## 兔抗 DEK 多克隆抗体

中文名称：兔抗 DEK 多克隆抗体

英文名称：Anti-DEK rabbit polyclonal antibody

别名：D6S231E

相关类别：一抗

抗原：DEK

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein with one SAP domain. This protein binds to cruciform and superhelical DNA and induces positive supercoils into closed circular DNA, and is also involved in splice site selection during mRNA processing. Chromosomal aberrations involving this region, increased expression of this gene, and the presence of antibodies against this protein are all associated with various diseases. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DEK

<b>Immunogen:</b>	Synthetic peptide of human DEK
<b>Full name:</b>	DEK proto-oncogene
<b>Synonyms :</b>	D6S231E
<b>SwissProt:</b>	P35659
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human gastric cancer
<b>IHC Recommend dilution:</b>	25-100

