

## 兔抗 DHRS7 多克隆抗体

中文名称：兔抗 DHRS7 多克隆抗体

英文名称： Anti-DHRS7 rabbit polyclonal antibody

别名： dehydrogenase/reductase 7; CGI-86; SDR34C1; retDSR4; retSDR4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DHRS7

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family, which has over 46,000 members. Members in this family are enzymes that metabolize many different compounds, such as steroid hormones, prostaglandins, retinoids, lipids and xenobiotics.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	DHRS7
<b>Immunogen:</b>	Fusion protein of human DHRS7
<b>Full name:</b>	dehydrogenase/reductase 7
<b>Synonyms:</b>	CGI-86; SDR34C1; retDSR4; retSDR4
<b>SwissProt:</b>	Q9Y394
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	38 kDa
<b>WB Positive control:</b>	Human kidney tissue lysate
<b>WB Recommended dilution:</b>	500-2000



