

兔抗 DIMT1 多克隆抗体

中文名称：兔抗 DIMT1 多克隆抗体

英文名称： Anti-DIMT1 rabbit polyclonal antibody

别名： DIM1 dimethyladenosine transferase 1 homolog; DIM1; DIMT1L; HUSSY5; HSA9761

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DIMT1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Specifically dimethylates two adjacent adenosines in the loop of a conserved hairpin near the 3'-end of 18 S rRNA in the 40S particle.
Applications:	ELISA, WB, IHC
Name of antibody:	DIMT1
Immunogen:	Full length fusion protein
Full name:	DIM1 dimethyladenosine transferase 1 homolog
Synonyms:	DIM1; DIMT1L; HUSSY5; HSA9761
SwissProt:	Q9UNQ2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	30-150

WB Predicted band size:	35 kDa
WB Positive control:	231 and Raji cell lysates
WB Recommended dilution:	500-2000



