

兔抗 DIPK2B 多克隆抗体

中文名称: 兔抗 DIPK2B 多克隆抗体

英文名称: Anti-DIPK2B rabbit polyclonal antibody

别名: divergent protein kinase domain 2B; DIA1R; CXorf36

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: DIPK2B

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

| | |
|--------------------------------|---|
| Background: | The function of this protein remains unknown. |
| Applications: | ELISA, IHC |
| Name of antibody: | DIPK2B |
| Immunogen: | Synthetic peptide of human DIPK2B |
| Full name: | divergent protein kinase domain 2B |
| Synonyms: | DIA1R; CXorf36; PRO3743; EPQL1862; bA435K1.1; 493057 8C19Rik |
| SwissProt: | Q9H7Y0 |
| IHC positive control: | Human liver cancer |
| IHC Recommend dilution: | 25-100 |

