

## 兔抗 DLL1 多克隆抗体

中文名称：兔抗 DLL1 多克隆抗体

英文名称： Anti-DLL1 rabbit polyclonal antibody

别名： ODF; OPGL; sOdf; CD254; OPTB2; RANKL; TNLG6B; TRANCE; hRANKL2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DLL1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	DLL1 is a human homolog of the Notch Delta ligand and is a member of the delta/serrate/jagged family. It plays a role in mediating cell fate decisions during hematopoiesis . It may play a role in cell-to-cell communication.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	DLL1
<b>Immunogen:</b>	Synthetic peptide of human DLL1
<b>Full name:</b>	delta like canonical Notch ligand 1
<b>Synonyms:</b>	DL1; Delta; DELTA1
<b>SwissProt:</b>	O00548
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human esophagus cancer and Human brain
<b>IHC Recommend dilution:</b>	10-50
<b>WB Predicted band size:</b>	78 kDa
<b>WB Positive control:</b>	Rat heart tissue lysate
<b>WB Recommended dilution:</b>	200-1000



