

兔抗 DNPEP 多克隆抗体

中文名称：兔抗 DNPEP 多克隆抗体

英文名称：Anti-DNPEP rabbit polyclonal antibody

别名：aspartyl aminopeptidase; DAP; ASPEP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DNPEP

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is an aminopeptidase which prefers acidic amino acids, and specifically favors aspartic acid over glutamic acid. It is thought to be a cytosolic protein involved in general metabolism of intracellular proteins. Several transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	DNPEP
Immunogen:	Fusion protein of human DNPEP
Full name:	aspartyl aminopeptidase
Synonyms:	DAP; ASPEP
SwissProt:	Q9ULA0

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human cervical cancer
IHC Recommend dilution:	30-150



