

兔抗 DPYSL3 多克隆抗体

中文名称：兔抗 DPYSL3 多克隆抗体

英文名称： Anti-DPYSL3 rabbit polyclonal antibody

别名： dihydropyrimidinase like 3; DRP3; ULIP; CRMP4; DRP-3; LCRMP; CRMP-4; ULIP-1

相关类别： 一抗

抗原： DPYSL3

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton. Plays a role in axon guidance, neuronal growth cone collapse and cell migration (By similarity).
Applications:	ELISA, IHC
Name of antibody:	DPYSL3
Immunogen:	Fusion protein of human DPYSL3
Full name:	dihydropyrimidinase like 3
Synonyms:	DRP3; ULIP; CRMP4; DRP-3; LCRMP; CRMP-4; ULIP-1
SwissProt:	Q14195
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer and Human esophagus cancer

IHC Recommend dilution: 50-100



