

## 兔抗 DRD3 多克隆抗体

中文名称：兔抗 DRD3 多克隆抗体

英文名称： Anti-DRD3 rabbit polyclonal antibody

别名： D3DR; ETM1; FET1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DRD3

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

This gene encodes the D3 subtype of the five (D1-D5) dopamine receptors. The activity of the D3 subtype receptor is mediated by G proteins which inhibit adenylyl cyclase. This receptor is localized to the limbic areas of the brain, which are associated with cognitive, emotional, and endocrine functions. Genetic variation in this gene may be associated with susceptibility to hereditary essential tremor 1. Alternative splicing of this gene results in transcript variants encoding different isoforms, although some variants may be subject to nonsense-mediated decay (NMD).

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DRD3
<b>Immunogen:</b>	Synthetic peptide of human DRD3
<b>Full name:</b>	dopamine receptor D3
<b>Synonyms :</b>	D3DR; ETM1; FET1
<b>SwissProt:</b>	P35462
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	30-150

