

兔抗 DUSP16 多克隆抗体

中文名称：兔抗 DUSP16 多克隆抗体

英文名称：Anti-DUSP16 rabbit polyclonal antibody

别名：MKP7; MKP-7

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DUSP16

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a mitogen-activated protein kinase phosphatase that is a member of the dual specificity protein phosphatase subfamily. These phosphatases inactivate their target kinases by dephosphorylating both the phosphoserine/threonine and phosphotyrosine residues. The encoded protein specifically regulates the c-Jun amino-terminal kinase (JNK) and extracellular signal-regulated kinase (ERK) pathways.
Applications:	ELISA, IHC
Name of antibody:	DUSP16
Immunogen:	Synthetic peptide of human DUSP16
Full name:	dual specificity phosphatase 16

Synonyms :	MKP7; MKP-7
SwissProt:	Q9BY84
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and human cervical cancer
IHC Recommend dilution:	25-100



