

## 兔抗 DUSP3 多克隆抗体

中文名称：兔抗 DUSP3 多克隆抗体

英文名称：Anti-DUSP3 rabbit polyclonal antibody

别名：VHR

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DUSP3

反应种属：Human, Mouse


标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the dual specificity protein phosphatase subfamily. This gene maps in a region that contains the BRCA1 locus which confers susceptibility to breast and ovarian cancer. Although DUSP3 is expressed in both breast and ovarian tissues, mutation screening in breast cancer pedigrees and in sporadic tumors was negative, leading to the conclusion that this gene is not BRCA1.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	DUSP3
<b>Immunogen:</b>	Synthetic peptide of human DUSP3
<b>Full name:</b>	dual specificity phosphatase 3

<b>Synonyms :</b>	VHR
<b>SwissProt:</b>	P51452
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer and human brain
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	20 kDa
<b>WB Positive control:</b>	Human fetal liver tissue
<b>WB Recommended dilution:</b>	200-1000



**mlbio** 梅联生物  
Good elisakit producers

