

兔抗 DYNLL2 多克隆抗体

中文名称：兔抗 DYNLL2 多克隆抗体

英文名称：Anti-DYNLL2 rabbit polyclonal antibody

别名：Dlc2; DNCL1B; RSPH22

抗原：DYNLL2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules. May play a role in changing or maintaining the spatial distribution of cytoskeletal structures (By similarity).
Applications:	ELISA, IHC
Name of antibody:	DYNLL2
Immunogen:	Full length fusion protein
Full name:	dynein, light chain, LC8-type 2
Synonyms :	Dlc2; DNCL1B; RSPH22

SwissProt:	Q96FJ2
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100

