

兔抗 E2F3 多克隆抗体

中文名称：兔抗 E2F3 多克隆抗体

英文名称：Anti-E2F3 rabbit polyclonal antibody

别名：E2F-3

相关类别：一抗

抗原：E2F3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of a small family of transcription factors that function through binding of DP interaction partner proteins. The encoded protein recognizes a specific sequence motif in DNA and interacts directly with the retinoblastoma protein (pRB) to regulate the expression of genes involved in the cell cycle. Altered copy number and activity of this gene have been observed in a number of human cancers. There are pseudogenes for this gene on chromosomes 2 and 17. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC

Name of antibody:	E2F3
Immunogen:	Synthetic peptide of human E2F3
Full name:	E2F transcription factor 3
Synonyms:	E2F-3
SwissProt:	O00716
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	25-100

