

兔抗 EEF2 多克隆抗体

中文名称：兔抗 EEF2 多克隆抗体

英文名称： Anti-EEF2 rabbit polyclonal antibody

别名： EF2; EF-2; EEF-2; SCA26

相关类别： 一抗

抗原： EEF2

储存： 冷冻 (-20℃)

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the GTP-binding translation elongation factor family. This protein is an essential factor for protein synthesis. It promotes the GTP-dependent translocation of the nascent protein chain from the A-site to the P-site of the ribosome. This protein is completely inactivated by EF-2 kinase phosphorylation.
Applications:	ELISA, WB, IHC
Name of antibody:	EEF2
Immunogen:	Fusion protein of human EEF2
Full name:	eukaryotic translation elongation factor 2
Synonyms:	EF2; EF-2; EEF-2; SCA26
SwissProt:	P13639

IHC positive control:	Human gastric cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	95 kDa
WB Positive control:	LoVo, HepG2 and MCF-7 cells
WB Recommended dilution:	500-2000



