

兔抗 CSPG4 多克隆抗体

- 中文名称：兔抗 CSPG4 多克隆抗体
- 英文名称：Anti-CSPG4 rabbit polyclonal antibody
- 别名：NG2; MCSP; MCSPG; MSK16; HMW-MAA; MEL-CSPG
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：CSPG4
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	A human melanoma-associated chondroitin sulfate proteoglycan plays a role in stabilizing cell-substratum interactions during early events of melanoma cell spreading on endothelial basement membranes. CSPG4 represents an integral membrane chondroitin sulfate proteoglycan expressed by human malignant melanoma cells.
Applications:	ELISA, WB, IHC
Name of antibody:	CSPG4
Immunogen:	Synthetic peptide of human CSPG4
Full name:	chondroitin sulfate proteoglycan 4
Synonyms :	NG2; MCSP; MCSPG; MSK16; HMW-MAA; MEL-CSPG
SwissProt:	Q6UVK1

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain
IHC Recommend dilution:	20-100
WB Predicted band size:	251 kDa
WB Positive control:	Mouse skin tissue
WB Recommended dilution:	500-2000



