

HL-1 完全培养基

基本信息

细胞名称	HL-1 完全培养基
细胞货号	ml-CC3026
细胞品牌	酶联生物
细胞规格	500ml
细胞描述	<p>HL-1 Cardiac Muscle Cell Line is an immortalized mouse cardiomyocyte cell line able to continuously divide and spontaneously contract while maintaining a differentiated cardiac phenotype. HL-1 can be serially passaged without losing its differentiated cardiac myocyte phenotype, including morphological, biochemical, and electrophysiological properties. HL-1 is used in a variety of model systems to address questions of cardiac biology at the cellular & molecular levels.</p>
细胞形态	液体
培养基成分	MEM+10%FBS+1%P/S
支原体检测	阴性
细胞生长	细胞生长良好, 形态正常
细胞货期	现货, 1周左右
储存条件	2~8°C, 避光储存

有效期	3个月
注意事项	使用时应注意无菌操作，避免污染。为保持本产品的使用效果，不宜长时间放置于室温或较高的温度环境中。冻融后，可能会有少量絮状物析出，不影响正常使用，超出保质期，必须放弃使用。

售后服务

细胞予重发

1. 细胞运输中遭遇的各种问题，细胞丢失瓶身破损、培养液严重漏液等，重发。
2. 收到细胞未开封，如出现污染状况，重发。
3. 收到细胞3天内，发现污染问题，经核实后，重发。
4. 常温发货的细胞静置2小时后，干冰冻存发货的细胞复苏2天后，绝大多数细胞未存活，经核实后，重发。
5. 常温发货的细胞静置22小时并且未开封或干冰冻存发货的细胞复苏2天后，出现污染，经核实后，重发。
6. 细胞活性问题，请在收到产品3天内给我们提出真实的实验结果，用台盼蓝染色法鉴定细胞活力，经核实后，重发。

细胞不重发

1. 客户操作造成细胞污染，不重发。
2. 客户严重操作失误致细胞状态不好，不重发。
3. 非我们推荐细胞培养体系致的细胞状态不好，不重发。
4. 细胞状态不好，未提供真实清晰的培养前3天的细胞状态照片，不重发。
5. 细胞培养时经其它处理导致细胞出现问题的，不重发。

6. 收到细胞发现问题与客服人员沟通的时间证明大于 3 天的, 不重发。

特别说明

客户买细胞就找**上海酶联生物**, 稳定传代, 无污染, 包存活, 提供整体课题外包服务, 光学成像, 流式实验, 电镜实验, 动物实验, 病理实验, 分子生物学实验, 细胞实验等, 严格把控产品质量, 所有细胞产品均有细胞鉴别、无菌检查、支原体检查, 为科研人员提供可靠放心的产品。